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ANALYSIS	
TYPE	RATE
Main W.Length 1	340
Main W.Length 2	450
METHOD	I

CORR.	
Y =	SLOPE
	1
X +	INTER
	0

CALIBRATION			
TYPE	Factor		
FACTOR -1750**			

NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum	6	40	6	32
Urine				
Plasma				
CSF				
Dialysis				
Other				

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Diluent ☐ 99:Dil1 ☒ 100:Dil2

PROZONE CHECK

	START	END	LIMIT(%)	
FIRST				
SECOND				☐ Low ☐ High
THIRD				☐ Low ☐ High

MONITOR

0 LEVEL POINT	1
SPAN	3

Page : 2

Item Name		#	ALAT
-----------	--	---	------

Auto Rerun SW

☐ On ☒ Off

Auto Rerun Condition (Absorbance)

Absorbance Range

Lower	☐ On ☒ Off
Higher	☐ On ☒ Off

Prozone Range ☐ On ☒ Off

Auto Rerun Range (Result)

	☐ On ☒ Off	☐ On ☒ Off
	Lower	Higher
Serum		450
Urine		
Plasma		
CSF		
Dialysis		
Other		

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NORMAL RANGE

	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum	#	#	#	#
Urine				
Plasma				
CSF				
Dialysis				
Other				

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Diluent ☐ 99:Dil1 ☒ 100:Dil2

PROZONE CHECK

	START	END	LIMIT(%)	
FIRST				
SECOND				☒ Low ☐ High
THIRD				☒ Low ☐ High

Page : 2

Item Name		#	ALB
-----------	--	---	-----

Auto Rerun SW

☐ On ☒ Off

Auto Rerun Range (Result)

	☐ On ☒ Off	☐ On ☐ Off
	Lower	Higher
Serum		6.0
Urine		
Plasma		
CSF		
Dialysis		
Other		

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Diluent ☐ 99:Dil1 ☒ 100:Dil2

PROZONE CHECK

	START	END	LIMIT(%)	
FIRST				☐ Low ☐ High ☐ Low ☐ High
SECOND				
THIRD				

MONITOR

0 LEVEL POINT	1
SPAN	3

Page : 2

Item Name	3	AST		
-----------	---	-----	--	--

Auto Rerun SW

☐ On ☒ Off

Auto Rerun Condition (Absorbance)

Absorbance Range

Lower	☐ On	☒ Off
Higher	☐ On	☒ Off

Prozone Range ☐ On ☒ Off

Auto Rerun Range (Result)

	☐ On ☒ Off	☒ On ☐ Off
	Lower	Higher
Serum		450
Urine		
Plasma		
CSF		
Dialysis		
Other		

Page : 3

*input calibrator value
entered by the user

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Diluent ☐ 99:Dil1 ☒ 100:Dil2

PROZONE CHECK

	START	END	LIMIT(%)	
FIRST				
SECOND				☐ Low ☐ High
THIRD				☐ Low ☐ High

MONITOR

0 LEVEL POINT	1
SPAN	3

Page : 2

Item Name	#	DBIL
-----------	---	------

Auto Rerun SW

☐ On ☒ Off

Auto Rerun Condition (Absorbance)

Absorbance Range

Lower	☐ On ☒ Off
Higher	☐ On ☒ Off

Prozone Range ☐ On ☒ Off

Auto Rerun Range (Result)

	☐ On ☒ Off	☐ On ☒ Off
	Lower	Higher
Serum		10.0
Urine		
Plasma		
CSF		
Dialysis		
Other		

Page : 3

*input calibrator value
entered by the user

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NORMAL RANGE

	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum	0.50	6.20	0.50	6.20
Urine				
Plasma				
CSF				
Dialysis				
Other				

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☐ 99:Dil1
 ☒ 100:Dil2

PROZONE CHECK

	START	END	LIMIT(%)	
FIRST				☐ Low ☐ High ☐ Low ☐ High
SECOND				
THIRD				

Page : 2

Item Name	#	BILT		
-----------	---	------	--	--

Auto Rerun SW

☒ On
 ☐ Off

Auto Rerun Range (Result)

	☐ On ☒ Off	☐ On ☐ Off
	Lower	Higher
Serum		150
Urine		
Plasma		
CSF		
Dialysis		
Other		

Auto Rerun Condition (Absorbance)

Absorbance Range

Lower	☐ On ☒ Off
Higher	☒ On ☐ Off

Prozone Range

☐ On ☒ Off

Page : 3

*input calibrator value
entered by the user

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CORR.

Y =	<table border="1"><tr><td>SLOPE</td></tr><tr><td>1</td></tr></table>	SLOPE	1	X +	<table border="1"><tr><td>INTER</td></tr><tr><td>0</td></tr></table>	INTER	0
SLOPE							
1							
INTER							
0							

CALIBRATION

TYPE	Factor
FACTOR 1158**	

NORMAL RANGE

	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum	10	90	10	66
Urine				
Plasma				
CSF				
Dialysis				
Other				

*input calibrator value
entered by the user

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Diluent ☐ 99:Dil1 ☒ 100:Dil2

PROZONE CHECK

	START	END	LIMIT(%)	
FIRST				☐ Low ☐ High ☐ Low ☐ High
SECOND				
THIRD				

MONITOR

0 LEVEL POINT	1
SPAN	3

Page : 2

Item Name	#	CA		
-----------	---	----	--	--

Auto Rerun SW

☒ On ☐ Off

Auto Rerun Condition (Absorbance)

Absorbance Range

Lower	☐ On	☒ Off
Higher	☒ On	☐ Off

Prozone Range ☐ On ☒ Off

Auto Rerun Range (Result)

	☐ On ☒ Off	☐ On ☐ Off
	Lower	Higher
Serum		230
Urine		
Plasma		
CSF		
Dialysis		
Other		

Page : 3

*input calibrator value
entered by the user

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NORMAL RANGE

	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum	3.88	6.10	3.88	6.10
Urine				
Plasma				
CSF				
Dialysis				
Other				

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LINEAR CHECK(%) ☐ 0

PROZONE CHECK

	START	END	LIMIT(%)	
FIRST				
SECOND				☒ Low ☐ High
THIRD				☒ Low ☐ High

Page : 2

Item Name	#	GLUC		
-----------	---	------	--	--

Auto Rerun SW

☒ On ☐ Off

Auto Rerun Range (Result)

	☒ On ☐ Off	☒ On ☐ Off
	Lower	Higher
Serum		500.0
Urine		
Plasma		
CSF		
Dialysis		
Other		

Auto Rerun Condition (Absorbance)

Absorbance Range

Lower	☐ On ☒ Off
Higher	☒ On ☐ Off

Prozone Range ☐ On ☒ Off

Page : 3

*input calibrator value
entered by the user

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Diluent ☐ 99:Dil1 ☒ 100:Dil2

PROZONE CHECK

	START	END	LIMIT(%)	
FIRST				
SECOND				☐ Low ☐ High
THIRD				☐ Low ☐ High

Page : 2

Item Name	#	GLUC
-----------	---	------

Auto Rerun SW

☐ On ☒ Off

Auto Rerun Range (Result)

	☐ On ☒ Off	☐ On ☒ Off
	Lower	Higher
Serum		20
Urine		
Plasma		
CSF		
Dialysis		
Other		

Auto Rerun Condition (Absorbance)

Absorbance Range

Lower	☐ On ☒ Off	
Higher	☐ On ☒ Off	

Prozone Range ☐ On ☒ Off

Page : 3

*input calibrator value
entered by the user

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CORR.					
Y =	<table border="1"><tr><td>SLOPE</td><td>INTER</td></tr><tr><td>1</td><td>0</td></tr></table> X +	SLOPE	INTER	1	0
SLOPE	INTER				
1	0				

CALIBRATION				
TYPE		Linear		
STANDARD				
#1	*	#4		
#2		#5		
#3		#6		

NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum	11.6	31.3	7.16	26.85
Urine				
Plasma				
CSF				
Dialysis				
Other				

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Prozone Range ☐ On ☒ Off

Auto Rerun Range (Result)

	☐ On ☒ Off	☒ On ☐ Off
	Lower	Higher
Serum		1000
Urine		
Plasma		
CSF		
Dialysis		
Other		

Page : 3

*input calibrator value
entered by the user

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ANALYSIS

TYPE	END
Main W.Length 1	570
Main W.Length 2	660
METHOD	O-Cresolphth

CORR.

Y =	<table border="1"><tr><td>SLOPE</td></tr><tr><td>1</td></tr></table>	SLOPE	1	X +	<table border="1"><tr><td>INTER</td></tr><tr><td>0</td></tr></table>	INTER	0
SLOPE							
1							
INTER							
0							

CALIBRATION

TYPE	Linear		
STANDARD			
#1	2,50	#4	
#2		#5	
#3		#6	

NORMAL RANGE

	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum	2,15	2,65	2,15	2,65
Urine				
Plasma				
CSF				
Dialysis				
Other				

Page : 1

Item Name	#	CA		
-----------	---	----	--	--

ASPIRATION

KIND ☐ Single ☐ Double

SAMPLE	VOLUME	
REAGENT1 VOL	3	μl
REAGENT2	300	

DATA PROCESS

READ

	START	END
MAIN	52	54
SUB		

ABSORBANCE LIMIT

LOW	0.000
HIGH	3.000

FACTOR

Blk correction	1
----------------	---

ENDPOINT LIMIT 2.000

LINEAR CHECK(%) 0

Dilution

Diluent ☐ 99:Dil1 ☐ 100:Dil2

PROZONE CHECK

	START	END	LIMIT(%)	
FIRST				☐ Low ☐ High ☐ Low ☐ High
SECOND				
THIRD				

MONITOR

0 LEVEL POINT	1
SPAN	3

Page : 2

Item Name	#	CA		
-----------	---	----	--	--

Auto Rerun SW

☐ On ☐ Off

Auto Rerun Condition (Absorbance)

Absorbance Range

Lower	☐ On	☐ Off
Higher	☐ On	☐ Off

Prozone Range ☐ On ☐ Off

Auto Rerun Range (Result)

	☐ On ☐ Off	☐ On ☐ Off
	Lower	Higher
Serum		5,00
Urine		
Plasma		
CSF		
Dialysis		
Other		

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NORMAL RANGE

	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum	44	100	44	88
Urine				
Plasma				
CSF				
Dialysis				
Other				

Page : 1

*input calibrator value

entered by the user

Item Name	#	Crea
-----------	---	------

ASPIRATION

KIND ☒ Single ☐ Double

SAMPLE	VOLUME	
REAGENT1 VOL	25	μl
REAGENT2	250	

Third Mix. ☐ OFF ☐ ON

R1 Blank ☒ Water Blk ☐ R1 Blank1

MONITOR

0 LEVEL POINT	1
SPAN	3

DATA PROCESS

READ

	START	END
MAIN	14	17
SUB		

ABSORBANCE LIMIT

LOW	0.000
HIGH	3.000

FACTOR

Blk correction	1
----------------	---

ENDPOINT LIMIT

LINEAR CHECK(%)

Dilution

Diluent ☐ 99:Dil1 ☒ 100:Dil2

PROZONE CHECK

	START	END	LIMIT(%)	
FIRST				
SECOND				☐ Low ☐ High
THIRD				☐ Low ☐ High

Page : 2

Item Name	#	CKMB
-----------	---	------

Auto Rerun SW

☒ On ☐ Off

Auto Rerun Range (Result)

	☐ On ☐ Off	☐ On ☐ Off
	Lower	Higher
Serum		900
Urine		
Plasma		
CSF		
Dialysis		
Other		

Auto Rerun Condition (Absorbance)

Absorbance Range

Lower	☐ On	☐ Off
Higher	☐ On	☐ Off

Prozone Range ☐ On ☐ Off

Page : 3

*input calibrator value
entered by the user

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ANALYSIS	
TYPE	END
Main W.Length 1	546
Main W.Length 2	-
METHOD	

CORR.					
Y =	<table border="1"><tr><td>SLOPE</td></tr><tr><td>1</td></tr></table> X + <table border="1"><tr><td>INTER</td></tr><tr><td>0</td></tr></table>	SLOPE	1	INTER	0
SLOPE					
1					
INTER					
0					

CALIBRATION			
TYPE		Linear	
STANDARD			
#1	357	#4	
#2		#5	
#3		#6	

NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum	202	414	142	389
Urine				
Plasma				
CSF				
Dialysis				
Other				

Page : 1

Item Name	#	UA		
-----------	---	----	--	--

ASPIRATION

KIND ☒ Single ☐ Double

SAMPLE	VOLUME	
REAGENT1 VOL	5	μl
REAGENT2	200	

Third Mix. ☐ OFF ☐ ON

R1 Blank ☐ Water Blk ☐ R1 Blank1

MONITOR

0 LEVEL POINT	1
SPAN	3

DATA PROCESS

READ

	START	END
MAIN	52	54
SUB		

ABSORBANCE LIMIT

LOW	0.000
HIGH	3.000

FACTOR

Blk correction	0
----------------	---

ENDPOINT LIMIT

LINEAR CHECK(%)

Dilution

Diluent ☒ 99:Dil1 ☐ 100:Dil2

PROZONE CHECK

	START	END	LIMIT(%)	
FIRST				
SECOND				☐ Low ☐ High
THIRD				☐ Low ☐ High

Page : 2

Item Name	#	UA		
-----------	---	----	--	--

Auto Rerun SW

☐ On ☐ Off

Auto Rerun Range (Result)

	☐ On ☐ Off	☐ On ☐ Off
	Lower	Higher
Serum		
Urine		
Plasma		
CSF		
Dialysis		
Other		

Auto Rerun Condition (Absorbance)

Absorbance Range

Lower	☐ On ☐ Off	☐ On ☐ Off
Higher	☐ On ☐ Off	☐ On ☐ Off

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CORR.					
Y =	<table border="1"><tr><td>SLOPE</td><td>INTER</td></tr><tr><td>1</td><td>0</td></tr></table> X +	SLOPE	INTER	1	0
SLOPE	INTER				
1	0				

CALIBRATION			
TYPE		Linear	
STANDARD			
#1	8.33	#4	
#2		#5	
#3		#6	

NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum	2.50	8.32	2.50	8.32
Urine				
Plasma				
CSF				
Dialysis				
Other				

Page : 1

*input calibrator value
entered by the user

Item Name	#	Urea
-----------	---	------

ASPIRATION

KIND ☒ Single ☐ Double

SAMPLE	VOLUME	
REAGENT1 VOL	3	μl
REAGENT2	300	

Third Mix. ☐ OFF ☐ ON

R1 Blank ☐ Water Blk ☒ R1 Blank1

MONITOR

0 LEVEL POINT	1
SPAN	3

DATA PROCESS

READ

	START	END
MAIN	26	30
SUB		

ABSORBANCE LIMIT

LOW	0.000
HIGH	3.000

FACTOR

Blk correction	1
----------------	---

ENDPOINT LIMIT

LINEAR CHECK(%)

Dilution

Diluent ☒ 99:Dil1 ☐ 100:Dil2

PROZONE CHECK

	START	END	LIMIT(%)	
FIRST				
SECOND				☐ Low ☐ High
THIRD				☐ Low ☐ High

Page : 2

Item Name	#	Urea
-----------	---	------

Auto Rerun SW

☒ On ☐ Off

Auto Rerun Range (Result)

	☐ On ☐ Off	☐ On ☐ Off
	Lower	Higher
Serum		35
Urine		
Plasma		
CSF		
Dialysis		
Other		

Auto Rerun Condition (Absorbance)

Absorbance Range

Lower	☐ On	☐ Off
Higher	☒ On	☐ Off

Prozone Range ☐ On ☒ Off

Page : 3

*input calibrator value
entered by the user

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CORR.

Y =	<table border="1"><tr><td>SLOPE</td></tr><tr><td>1</td></tr></table>	SLOPE	1	X +	<table border="1"><tr><td>INTER</td></tr><tr><td>0</td></tr></table>	INTER	0
SLOPE							
1							
INTER							
0							

CALIBRATION

TYPE	Linear		
<u>STANDARD</u>			
#1	60	#4	
#2		#5	
#3		#6	

NORMAL RANGE

	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum	65	85	65	85
Urine				
Plasma				
CSF				
Dialysis				
Other				

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Diluent ☐ 99:Dil1 ☒ 100:Dil2

PROZONE CHECK

	START	END	LIMIT(%)	
FIRST				
SECOND				☐ Low ☐ High
THIRD				☐ Low ☐ High

MONITOR

0 LEVEL POINT	1
SPAN	3

Page : 2

Item Name	4	TP		
-----------	---	----	--	--

Auto Rerun SW

☐ On ☒ Off

Auto Rerun Condition (Absorbance)

Absorbance Range

Lower	☐ On	☒ Off
Higher	☐ On	☒ Off

Prozone Range ☐ On ☒ Off

Auto Rerun Range (Result)

	☐ On ☒ Off	☐ On ☒ Off
	Lower	Higher
Serum		850
Urine		
Plasma		
CSF		
Dialysis		
Other		

Page : 3

*input calibrator value
entered by the user

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NORMAL RANGE

	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum		1,92		1,92
Urine				
Plasma				
CSF				
Dialysis				
Other				

Page : 1

Item Name	#	TRIG		
-----------	---	------	--	--

ASPIRATION

KIND ☐ Single ☐ Double

SAMPLE	VOLUME	
REAGENT1 VOL	3	μl
REAGENT2	300	

DATA PROCESS

READ

	START	END
MAIN	52	54
SUB		

ABSORBANCE LIMIT

LOW	0.000
HIGH	3.000

FACTOR

Blk correction	1
----------------	---

ENDPOINT LIMIT 2.000

LINEAR CHECK(%) 0

Dilution

Diluent ☐ 99:Dil1 ☐ 100:Dil2

PROZONE CHECK

	START	END	LIMIT(%)	
FIRST				☐ Low ☐ High ☐ Low ☐ High
SECOND				
THIRD				

MONITOR

0 LEVEL POINT	1
SPAN	3

Page : 2

Item Name	#	TRIG		
-----------	---	------	--	--

Auto Rerun SW

☐ On ☐ Off

Auto Rerun Condition (Absorbance)

Absorbance Range

Lower	☐ On	☐ Off
Higher	☐ On	☐ Off

Prozone Range ☐ On ☐ Off

Auto Rerun Range (Result)

	☐ On ☐ Off	☐ On ☐ Off
	Lower	Higher
Serum		
Urine		
Plasma		
CSF		
Dialysis		
Other		

Page : 3

*input calibrator value
entered by the user

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ANALYSIS

TYPE	END
Main W.Length 1	505
Main W.Length 2	-
METHOD	Thiocyanate

CORR.

Y =	SLOPE	X +	INTER
	1		0

CALIBRATION

TYPE	Linear		
STANDARD			
#1	100	#4	
#2		#5	
#3		#6	

NORMAL RANGE

	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum	95	108	95	108
Urine				
Plasma				
CSF				
Dialysis				
Other				

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Auto Rerun Condition (Absorbance)

Absorbance Range		
Lower	☐ On ☐ Off	☐ On ☐ Off
Higher	☐ On ☐ Off	☐ On ☐ Off

Prozone Range	☐ On	☐ Off
---------------	--------------------------	---------------------------

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NORMAL RANGE

	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum	3.7	6.7	3.7	6.7
Urine				
Plasma				
CSF				
Dialysis				
Other				

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ANALYSIS	
TYPE	RETE
Main W.Length 1	405
Main W.Length 2	600
METHOD	IFCC

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NORMAL RANGE				
	MALE		FEMALE	
	LOW	HIGH	LOW	HIGH
Serum	40.2	129	34.8	104.4
Urine				
Plasma				
CSF				
Dialysis				
Other				

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Prozone Range ☐ On ☒ Off

Auto Rerun Range (Result)

	☐ On ☒ Off	☐ On ☒ Off
	Lower	Higher
Serum		
Urine		
Plasma		
CSF		
Dialysis		
Other		

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