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10 Feb 10
21:50:56
PAGE 1

CX4 PRO USER-DEFINED CHEMISTRIES

USER ID : 36
Chemistry Name: PALT
Test Name: PALT

Reaction Type: RATE 1
Reaction Direction: Decreasing
Units: U/L

Calculation Factor: 3599
Math Model: Linear
Cal Time Limit: 0 Hrs
No. of Calibrators: 0

Decimal Precision: X

Primary Wavelength: 340 nm Secondary Wavelength: 410 nm

Sample Volume: 25 µl CALIBRATORS MULTIPOINT SPAN

Primary Inject Rgt: -----

A : 250 µl
None : 0 µl

Secondary Inject Rgt: -----

None : 0 µl
Add Time: 0 sec

| | |
|---------------------|---------------------|
| REAGENT BLANK | REACTION |
| ----- | ----- |
| Start Read: 10 sec | Start Read: 60 sec |
| End Read: 30 sec | End Read: 180 sec |
| Low ABS Limit: -1.5 | Low ABS Limit: -1.5 |
| High ABS Limit: 1.5 | High ABS Limit: 1.5 |

| | |
|-----------------------|-----------------------|
| USABLE RANGE | SUBSTRATE DEPLETION |
| ----- | ----- |
| Lower Limit: 0.00 | Initial Rate: -99.999 |
| Upper Limit: 99999.00 | Delta ABS: 0.3 |

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USER ID : 50
Chemistry Name: PALB
Test Name: PALB

Calculation Factor: 0

Reaction Type: Endpoint
Reaction Direction: Increasing
Units: g/L

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Decimal Precision: X

Primary Wavelength: 600 nm Secondary Wavelength: 700 nm

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10 Feb 10
21:50:56
PAGE 1

CX4 PRO USER-DEFINED CHEMISTRIES

USER ID : 37
Chemistry Name: PAST
Test Name: PAST

Calculation Factor: 3142

Reaction Type: RATE 1
Reaction Direction: Decreasing
Units: U/L

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Decimal Precision: X

Primary Wavelength: 340 nm Secondary Wavelength: 410 nm

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10 Feb 10
21:51:23
PAGE 1

CX4 PRO USER-DEFINED CHEMISTRIES

USER ID : 51
Chemistry Name: PBIT
Test Name: PBIT
Calculation Factor: 600

Reaction Type: Endpoint
Reaction Direction: Increasing
Units: umol/L
Math Model: Linear
Cal Time Limit: 0 Hrs
No. of Calibrators: 0

Decimal Precision: X.X

Primary Wavelength: 560 nm
Secondary Wavelength: 700 nm

Sample Volume: 25 µl
Primary Inject Rgt: CALIBRATORS
A : 250 µl
None : 0 µl
Secondary Inject Rgt: MULTIPOINT SPAN
B : 25 µl
Add Time: 16 sec

| REAGENT BLANK | REACTION |
|---------------------|---------------------|
| Start Read: 250 sec | Start Read: 300 sec |
| End Read: 300 sec | End Read: 330 sec |
| Low ABS Limit: -1.5 | Low ABS Limit: -1.5 |
| High ABS Limit: 1.5 | High ABS Limit: 1.5 |

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Test Name: PBID

Calculation Factor: 600

Reaction Type: Endpoint
Reaction Direction: Increasing
Units: umol/L

Math Model: Linear
Cal Time Limit: 0 Hrs
No. of Calibrators: 0

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USER ID : 42
Chemistry Name: PGGT
Test Name: PGGT

Calculation Factor: 3076

Reaction Type: RATE 1
Reaction Direction: Increasing
Units: U/L

Math Model: Linear
Cal Time Limit: 0 Hrs
No. of Calibrators: 0

Decimal Precision: X

Primary Wavelength: 410 nm
Secondary Wavelength: 560 nm

Sample Volume: 25 μ l
Primary Inject Rgt:
A : 250 μ l
None : 0 μ l
Secondary Inject Rgt:
None : 0 μ l
Add Time: 0 sec

CALIBRATORS . MULTIPOINT SPAN

| REAGENT BLANK | REACTION |
|---------------------|---------------------|
| Start Read: 10 sec | Start Read: 60 sec |
| End Read: 30 sec | End Read: 180 sec |
| Low ABS Limit: -1.5 | Low ABS Limit: -1.5 |
| High ABS Limit: 1.5 | High ABS Limit: 1.5 |

| USABLE RANGE | SUBSTRATE DEPLETION |
|-----------------------|----------------------|
| Lower Limit: 0.00 | Initial Rate: 99.999 |
| Upper Limit: 99999.00 | Delta ABS: 1.5 |

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10 Feb 10
21:51:49
PAGE 1

CX4 PRO USER-DEFINED CHEMISTRIES

USER ID : 54
Chemistry Name: PGLG
Test Name: PGLG

Calculation Factor: 0

Reaction Type: Endpoint
Reaction Direction: Increasing
Units: mmol/L

Math Model: Linear
Cal Time Limit: 336 Hrs
No. of Calibrators: 1

Decimal Precision: X.X

Primary Wavelength: 340 nm Secondary Wavelength: 700 nm

Sample Volume: 3 µl
Primary Inject Rgt:
A : 300 µl
None : 0 µl
Secondary Inject Rgt:
None : 0 µl
Add Time: 0 sec

| CALIBRATORS | MULTIPOINT SPAN |
|-------------|-----------------|
| #1: 5.55 | |

REAGENT BLANK

Start Read: 250 sec
End Read: 300 sec
Low ABS Limit: -1.5
High ABS Limit: 1.5

USABLE RANGE

Lower Limit: 0.00
Upper Limit: 9999.99

REACTION

Start Read: 300 sec
End Read: 330 sec
Low ABS Limit: -1.5
High ABS Limit: 1.5

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Test Name: PGLU

Calculation Factor: 0

Reaction Type: Endpoint
Reaction Direction: Increasing
Units: nmol/L

Math Model: Linear
Cal Time Limit: 336 Hrs
No. of Calibrators: 1

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10 Feb 10
21:51:48
PAGE 1

CX4 PRO USER-DEFINED CHEMISTRIES

USER ID : 44
Chemistry Name: PAFE
Test Name: PAFE
Calculation Factor: 0
Reaction Type: Endpoint 2
Reaction Direction: Increasing
Units: umol/L
Math Model: Linear
Cal Time Limit: 72 Hrs
No. of Calibrators: 2
Decimal Precision: X.X

Primary Wavelength: 560 nm Secondary Wavelength: 600 nm

Sample Volume: 3 µl
Primary Inject Rgt:
A : 200 µl
None : 0 µl
Secondary Inject Rgt:
C : 6 µl
Add Time: 560 sec

| CALIBRATORS | MULTIPOINT SPAN |
|-------------|-----------------|
| #1: 0.00 | 1-2: 0.000 |
| #2: 35.70 | |

| REAGENT BLANK | REACTION |
|---------------------|---------------------|
| Start Read: 300 sec | Start Read: 600 sec |
| End Read: 512 sec | End Read: 700 sec |
| Low ABS Limit: -1.5 | Low ABS Limit: -1.5 |
| High ABS Limit: 1.5 | High ABS Limit: 1.5 |

| USABLE RANGE | SUBSTRATE DEPLETION |
|-----------------------|----------------------|
| Lower Limit: 0.00 | Initial Rate: 99.999 |
| Upper Limit: 99999.00 | Delta ABS: 1.5 |

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Test Name: PURE

Calculation Factor: 0

Reaction Type: RATE 1
Reaction Direction: Decreasing
Units: mmol/L

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10 Feb 10
21:51:23
PAGE 1

CX4 PRO USER-DEFINED CHEMISTRIES

USER ID : 46
Chemistry Name: PTP
Test Name: PTP

Calculation Factor: 0

Reaction Type: Endpoint 2
Reaction Direction: Increasing
Units: g/L

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Decimal Precision: X

Primary Wavelength: 520 nm Secondary Wavelength: 700 nm

Sample Volume: 6 µl CALIBRATORS MULTIPOINT SPAN

Primary Inject Rgt: -----

A : 300 µl #1: 60.00

None : 0 µl

Secondary Inject Rgt:

None : 0 µl

Add Time: 0 sec

| REAGENT BLANK | REACTION |
|---------------------|---------------------|
| Start Read: 250 sec | Start Read: 400 sec |
| End Read: 300 sec | End Read: 430 sec |
| Low ABS Limit: -1.5 | Low ABS Limit: -1.5 |
| High ABS Limit: 1.5 | High ABS Limit: 1.5 |

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11 Feb 10
14:50:40
PAGE 1

CX7 PRO USER-DEFINED CHEMISTRIES

USER ID : 18
Chemistry Name: PTRI
Test Name: PTRI
Calculation Factor: 0
Reaction Type: Endpoint
Reaction Direction: Decreasing
Units: mmol/L
Math Model: Linear
Cal Time Limit: 336 Hrs
No. of Calibrators: 1
Decimal Precision: X.XX
Primary Wavelength: 520 nm
Secondary Wavelength: 700 nm
Sample Volume: 3 µl
Primary Inject Rgt:
A : 300 µl
None : 0 µl
Secondary Inject Rgt:
None : 0 µl
Add Time: 0 sec

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Decimal Precision: X.XX

Primary Wavelength: 340 nm Secondary Wavelength: 560 nm

Sample Volume: 3 µl
Primary Inject Rgt:
A : 300 µl
None : 0 µl
Secondary Inject Rgt:
None : 0 µl
Add Time: 0 sec

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10 Feb 10
21:51:48
PAGE 1

CX4 PRO USER-DEFINED CHEMISTRIES

USER ID : 47
Chemistry Name: PCHO
Test Name: PCHO

Calculation Factor: 0

Reaction Type: Endpoint 2
Reaction Direction: Increasing
Units: nmol/L

Math Model: Linear
Cal Time Limit: 72 Hrs
No. of Calibrators: 1

Decimal Precision: X.X

Primary Wavelength: 520 nm
Secondary Wavelength: 700 nm

Sample Volume: 3 µl
Primary Inject Rgt:
A : 300 µl
None : 0 µl
Secondary Inject Rgt:
None : 0 µl
Add Time: 0 sec

CALIBRATORS
#1: 5.17

MULTIPOINT SPAN

| REAGENT BLANK | REACTION |
|---------------------|---------------------|
| Start Read: 256 sec | Start Read: 320 sec |
| End Read: 320 sec | End Read: 384 sec |
| Low ABS Limit: -1.5 | Low ABS Limit: -1.5 |
| High ABS Limit: 1.5 | High ABS Limit: 1.5 |

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10 Feb 10
21:50:56
PAGE 1

CX4 PRO USER-DEFINED CHEMISTRIES

USER ID : 43
Chemistry Name: PALP
Test Name: PALP
Calculation Factor: 4283

Reaction Type: RATE 1
Reaction Direction: Increasing
Units: U/L
Math Model: Linear
Cal Time Limit: 0 Hrs
No. of Calibrators: 0

Decimal Precision: X

Primary Wavelength: 410 nm Secondary Wavelength: 560 nm

Sample Volume: 6 µl
Primary Inject Rgt: CALIBRATORS MULTIPOINT SPAN
A : 300 µl
None : 0 µl
Secondary Inject Rgt: -----
None : 0 µl
Add Time: 0 sec

| | |
|------------------------|------------------------------|
| REAGENT BLANK ----- | REACTION ----- |
| Start Read: 10 sec | Start Read: 60 sec |
| End Read: 30 sec | End Read: 180 sec |
| Low ABS Limit: -1.5 | Low ABS Limit: -1.5 |
| High ABS Limit: 1.5 | High ABS Limit: 1.5 |
| USABLE RANGE ----- | SUBSTRATE DEPLETION ----- |
| Lower Limit: 0.00 | Initial Rate: 99.999 |
| Upper Limit: 99999.00 | Delta ABS: 1.5 |

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